

Affordable Wheelchair Lifts

Customer Pick-Up Instructions

Please read this if you are planning to come and pick up your lift at our manufacturing facility.

We are excited you will be acquiring your new Affordable Wheelchair Lift soon and look forward to greeting you upon your arrival at our facility to pick up your new lift. We hope it provides you and/or your loved ones with many years of improved mobility.

A few things to remember when preparing to pick up your lift at our facility:

- We are located in the: Business Development Centre
147 Mill Ridge Road
Suite 112
Lynchburg, VA. 24502
- Your primary contact person is: Cherrie Murphy
434-214-0204
- Your secondary contact person is: Richard Eagen
804-391-4002
- We have a forklift available to assist in loading the lift onto your vehicle or trailer.
- Please ensure you arrive during the scheduled pick up hours to ensure we have staff available to assist in loading your lift.

Transporting Your Lift:

Ensure your vehicle or trailer is an appropriate size to accommodate your lift and column.

Bring “tie downs” to secure your lift and shipping blankets to protect it. If it might rain you may consider bringing tarps and a means to secure them.

The standard lift platform is 52” x 43.5” x 40” and weighs around 200 lbs. If you have ordered a customized platform then these dimensions will be different. The length of the

lift's column varies from customer to customer. Depending on what options you have chosen there may be a few boxes with additional parts.

Many lifts are short enough to travel fully assembled. If possible, this is most convenient.

- The best way to transport your lift is by using a trailer with a ramp and a bed at least 60" wide. This will allow many lifts to use wheels to roll onto and off of the trailer and to arrive pre-assembled and ready to unwrap and plug in. Two men should be able to handle this.
- The next best way to transport your lift is in the back of a larger pickup truck *with no camper-shell*. The fully assembled lift can be laid on its back and placed on a pallet. Two or three men can unload it at the destination, set it upright and roll it into position.
- Another way to transport your lift is on a smaller trailer. The lift will not be as easy to roll on or off and may have to be disassembled to ship and reassembled at the destination.
- A final way to transport the lift is in a van or other vehicle. **Measure carefully to insure the lift will actually fit into the vehicle.** The lift will have to be disassembled to ship and reassembled at the destination.